
Granulation, Coating, Spheronization, and Mixing Equipment

ECX【pi:z】Ps for Lab

Equipped with spraying, temperature control, and dust collection functions



- Uniform coating achieved with a dedicated nozzle
- Suitable for temperature-sensitive products
- Prevents material adhesion inside the equipment through dust collection

Sealed cover designed for pharmaceutical applications



- Designed to prevent material leakage beyond the equipment
- Designed to prevent operator exposure to fine particles
- Designed to prevent external cross-contamination

ECX4Ps	
Overall Machine Dimensions	W825 × D950 × H1000 mm
Machine Weight	Approx. 400 kg
Charge Amount	Approx. 2L



ECX4Ps

Ideal for small-scale prototyping under temperature-controlled, contamination-free conditions

Tipton USA Inc.

Patent | JP4740419 JP5603808 JP7923007

Made in Japan

Head Office
3-19-21 Toyoda, Minami-Ku, Nagoya, Aichi Prefecture 457-8566 Japan
Tel : +81-52-692-6666
Fax : +81-52-692-9445

Headquarter-Sales
3-19-21 Toyoda, Minami-Ku, Nagoya, Aichi Prefecture 457-8566 Japan
Tel : +81-52-692-7175
Fax : +81-52-692-0249

Tobishima Factory
3-25-1 Odakara, Tobishima-Mura, Ama-Gun, Aichi Prefecture 490-1438 Japan
Tel : +81-567-56-7500
Fax : +81-567-56-7513

Tipton USA Inc.
2450 W Broadway Rd. Suite 108, Mesa, AZ 85202, USA
Tel : +1-480-795-2647

Quality Assurance Section for Product Quality
Tel : +81-567-56-7503
Fax : +81-567-56-7516

Engineering Section for Machine Maintenance
Tel : +81-567-56-7504
Fax : +81-567-56-7514

E-mail
sales-department@tipton.co.jp
sales-department@tiptonusa.com (North America)

Note) All unapproved reproduction of a report a photograph, etc are forbidden found in this catalog.

https://www.tipton.co.jp/

Catalog Download

YouTube

Access our social media for more!